

# Giraffe OmniBed



Giraffe\* OmniBed simplifies care with a full-featured incubator and a radiant warmer in one exceptional product.

## Overview:

### **Get the Best of Both Worlds—One Baby. One Bed.**

- A full-featured incubator and a radiant warmer in one. With the touch of a button, the OmniBed converts from a full-featured incubator to a radiant warmer. This helps reduce the stress to the baby and increase access for you.
- Integrated features create an all-in-one neonatal care station.
- In either mode, the Giraffe OmniBed combines advanced technology, exceptional Performance, and innovative design to provide an excellent, developmentally supportive, family-centered critical care environment for your high acuity newborns.

## Benefits:

### **Flexibility lets you respond quickly to the changing needs of patients.**

The Giraffe\* OmniBed is a combination system that brings together exceptional features of both a traditional incubator and a radiant warmer, so it truly can be at the center of the workspace for your premature patients.

### **We designed the workspace to make your job easier.**

With so much of what you need right at your fingertips, you will be able to provide the quality care you always have while dedicating even more time to your patients and their families.

## **Technology:**

### **The Giraffe\* OmniBed offers the following features:**

- The patented Baby Susan\* rotating mattress minimizes unnecessary baby stimulation. It rotates 360° and can slide out so you can easily position your patients for all types of procedures without disturbing them.
- The layered Giraffe pressure-diffusing mattress provides a soft support surface that can help relieve pressure points and preserve skin integrity.
- Large, drop down, removable doors and centralized control screen placement provide quick and easy patient access from either side of the bed.
- Elevating base offers a broad adjustable range to ensure a comfortable working height for procedures or family interaction while sitting down.
- Integrated scale allows quick and accurate assessment and trending of baby's growth and development.

## **Specifications:**

### **Dimensions**

- Maximum height: Canopy down: 70 in (178 cm) Canopy up: 94 in (239 cm)
- Minimum height: Canopy down: 58 in (147 cm) Canopy up: 82 in (208 cm)
- Footprint: 45 in x 26 in (114 cm x 66 cm)
- Weight: 284 lbs. (129 kg)
- Floor to mattress height: 32 in to 44 in (81 cm to 111 cm)
- Mattress size: 26 in x 19 in (66 cm x 48 cm)
- Access door size: 9 in x 31 in (23 cm x 78 cm)
- Drawer size: 19 in x 20 in (48 cm x 50 cm)
- Drawer depth: 8 in (20 cm)

### **Physical Characteristics**

- Mattress tilt angle: 12,° continuously variable
- Tubing access ports: 8
- Servo humidity control range: 30-95%
- relative humidity in 5% increments
- Alarm sound level: Adjustable audible levels

### **In-Bed Scale Performance**

- Accuracy:  $\pm 0.35$  oz (10 g)
- Range: 0.66 lbs to 17.6 lbs (300 g to 8 kg)