

# FlexiTrunk™ Midline Interface

- Foam block holder keeps the tubing in place with the infant's movement
- Capped pressure port for pressure monitor line connection
- Two models available with the standard model (BC18x) for connection to the F&P Bubble CPAP System and the universal model (BC19x) to connect to other circuits
- Patented glider technology for easy prong or mask repositioning and keeping straps away from the infant's face

Product Code	Description	Quantity
BC180-05	Nasal Tubing 50mm	5/box
BC181-05	Nasal Tubing 70mm	5/box
BC182-05	Nasal Tubing 100mm	5/box
BC190-05	Nasal Tubing Universal 50mm	5/box
BC191-05	Nasal Tubing Universal 70mm	5/box
BC192-05	Nasal Tubing Universal 100mm	5/box