

# Infant Nasal Prongs

**Product Code: See table below**

- Anatomically curved for comfortable fit
- Septum cut-away design to help prevent septum injury
- Soft, silicone, pliable and gentle on nares

| Part No.  | Description  | Quantity |
|-----------|--|----------|
| BC3020-10 | Nasal prong 3.0mm nare diameter/2.0mm septal width | 10/box   |
| BC3520-10 | Nasal prong 3.5mm nare diameter/2.0mm septal width | 10/box   |
| BC4030-10 | Nasal prong 4.0mm nare diameter/3.0mm septal width | 10/box   |
| BC4540-10 | Nasal prong 4.5mm nare diameter/4.0mm septal width | 10/box   |
| BC5040-10 | Nasal prong 5.0mm nare diameter/4.0mm septal width | 10/box   |
| BC5050-10 | Nasal prong 5.0mm nare diameter/5.0mm septal width | 10/box   |
| BC5550-10 | Nasal prong 5.5mm nare diameter/5.0mm septal width | 10/box   |
| BC5560-10 | Nasal prong 5.5mm nare diameter/6.0mm septal width | 10/box   |
| BC6060-10 | Nasal prong 6.0mm nare diameter/6.0mm septal width | 10/box   |
| BC6070-10 | Nasal prong 6.0mm nare diameter/7.0mm septal width | 10/box   |
| BC6570-10 | Nasal prong 6.5mm nare diameter/7.0mm septal width | 10/box   |